



EZ Way, Inc.

QUICK REFERENCE GUIDE eZ EFFORT



This quick Reference guide has been created to assist healthcare professionals and other users with only the functional aspects of the system. The guide does not replace the general safety warnings, cautions and recommendations provided in the product instructions for using the eZ Effort Lift®.



REFER TO THE FULL OPERATOR INSTRUCTIONS PRIOR TO USE.


Product Overview



Attach sling loops to the sling hanger hooks using the same length and color of loop on each side.




Use the hand control or secondary handle control to operate.





For emergency lowering only, turn the red collar on the actuator clockwise.

1. Fit the appropriate EZ Way sling to the person based on size and need. Confirm that the sling is properly positioned and that loop configurations match the transfer requirements.
2. Position the lift in front of the person, open the lift legs, if necessary, to get around obstacles. Ensure rear casters remain unlocked during transfer.
3. Lower the boom and position the nose cone over person's torso.
4. With the sling hangers parallel to the person's body, attach identical loop lengths and colors on both sides to the sling hangers.
5. Using the hand control or secondary control panel, raise the person to the needed height for the transfer.
6. Close the legs, if needed, and transfer the person to desired location.
7. Position the person over the target surface and lower to a sitting position or supine position, dependent on the transfer.
8. Lower the boom enough to remove the sling loops from the sling hanger hooks and move the unit away from the person.
9. Remove the sling from around the person, unless the sling is designed to remain under the person.



Leg Spreading - EL1000 Lift



Training Videos for
eZ Effort Lifts
<https://www.ezlifts.com/ezeffortliftqr>