

## **EZ REPOSITIONER**

# The EZ Repositioner allows one person to reposition a patient up or down in bed easily and pain-free.

With the EZ Repositioner, the patient's entire body is moved at the same time without being touched by the caregiver.

The EZ Repositioner easily attaches to most hospital beds using simple tools. A sheet is fed from a roller at the foot of the bed to a roller at the head of the bed. As the handle turns the roller, the sheet moves the patient toward the head of the bed without any lifting by the nurse or caregiver. By moving the handle to the roller at the foot of the bed, the repositioning sheet can be rewound for additional use.



Roller

Top Sheet

Roller



Bracket

Hospital Bed

**Bracket** 

Handle



Just **one caregiver** can easily reposition a patient with the turn of the crank. **5 yards of fabric** allow multiple repositionings. **Easy-to-clean** fabric resists stains.

### **Nurse & Caregiver Benefits**

- Virtually eliminates the risk of lower back injuries that result from repositioning the patient up in bed.
- Reduces costs and increases productivity by allowing a single nurse or caregiver to move the patient.

#### **Employer Benefits**

- Reduces lost work days due to back injuries associated with repositioning patients.
- Reduces workers' compensation injuries & insurance costs.

#### **Patient Benefits**

- Reduces shearing and friction to the patient's skin.
- Reduces the pain and discomfort associated with repositioning a patient in bed.
- Reduces the chance of injury to immobilized patients or patients with fragile or brittle bone conditions.
- Increases the patient's and caregiver's comfort and peace of mind by improving the quality of the care.

#### **Specifications**

All metal is steel. The sheet is an anti-bacterial, flame-retardant, anti-static, non-allergenic, tear-resistant material. EZ Way offers various brackets to fit a variety of bed makes and models.

#### Cleaning

Soils and stains are removed by wiping the fabric with neutral suds and lukewarm water. Rinse with water. Hard to clean spots will require the use of a standard housekeeping liquid cleaner and/or a soft bristle brush. Pre-soaking may be required. Do not use harsh cleaners or solvents.

Form #2-320 Rev. 08/08/2018