



EZ Way, Inc.

*"Your Partner in
Safe Patient Handling"*

VICTORIA

The RG9 Victoria is a freestanding tub that allows the caregiver 360° unobstructed access to the resident.



RG9D - Victoria Freestanding Bathing System

The low threshold and wide side-entry door allows easy access for residents with various degrees of mobility and the built-in contoured seat allows the resident to sit comfortably and safely in a slightly reclined position.

The sleek exterior design of the RG9 Victoria bathing system will give your bathing area an upscale look.

Caregivers will enjoy the RG9 Victoria for its ease of use and the improved quality of care they can provide the resident. Featuring a full range of standard and optional features, the RG9 Victoria bathing system is designed to enhance the bathing experience for the resident.



Easy to Use Control Panel with water temperature dial built in.



Shower Wand included as standard.



Secondary Locking Latch for added safety ideal for patients inclined to let themselves out of the tub.



The built-in disinfection system cabinet is conveniently located on the front side of the tub for easy access. The lockable cabinet provides a secure place for cleaning supplies and disinfectants.



The air jets on the Air Spa System feature spring loaded check valves for back flow prevention, which eliminates the extra step to disinfect the air lines.

"Your Partner in Safe Patient Handling"

Standard Features

- Thermostatic Mixing Valve with anti-scald protection that dispenses 15 gallons per minute for fast tub fill.
- Built-in Disinfection System for disinfection of the bath's external surfaces minimizes the risk of cross-infection.
- Handheld Shower Wand for hair and upper body washing, which can be used at the same time the tub is filling.
- Analog Thermometer for monitoring tub-fill and showerwand temperatures.
- Lap and Shoulder Harness for resident safety.
- Secondary Door Lock to prevent accidental door opening when in use.
- 3/4" Supply Lines constructed of stainless steel.
- 24-hour Service Department.

Why CleanRane™ Air Spa System?

The CleanRane™ Air Spa System enhances the therapeutic and relaxing effects of bathing and provides a pleasant experience for the resident. Hydromassage provides millions of preheated air bubbles in the water which increases stimulation of peripheral blood circulation, reduces pain and inflammation in the joints, gently massages away dirt, and induces deeper relaxation. The hygienic, low maintenance CleanRane™ Air Spa System features 12 air jets with spring loaded check valves for back flow prevention. Unlike most other whirlpool and air spa systems, the CleanRane™ Air Spa System does not allow water to enter the air lines, which can become a breeding ground for bacteria and mold. The motor is quiet, efficient, and has a ceramic heating element that blows warm air to help maintain the water temperature and warm the surface of the tub.

Technical Specifications

Length	60.75"
Width	33"
Height (to top of console)	43" (lowest setting)
Height (to top of door)	36.5" (lowest setting)
Adjustable Leveling Legs	2" max.
Seat Width (interior)	25"
Seat Height	18" (lowest setting)
Threshold Height	11" (lowest setting)
Transfer Device Access	5" H x 44" W
Weight Capacity (occupant)	600 lbs
Water Capacity (maximum)	52 gallons
Fill Time	3 minutes, 28 seconds
Drain Time	2 minutes, 30 seconds
Water Supply Inlet	3/4" diameter
Drain Size	1.5"
Electrical Supply	115V AC; 60Hz; 15 amp
Door Swing	Outward
Door Configuration	Right or Left Hinge



The wide side-entry door and unobstructed mobile lift access will allow the resident to enter and exit the bathtub safely and easily.

Warranty: 30 years for external finishes and composite shell, 1 year for all moving and stationary parts, and lifetime on the door seal.

Certifications / Standards: A112.19.7-2012 / CSA B45.10-12 / UL 1795 / CSA C22.2 no. 218.2 / IAPMO Z124-2011 / CSA B45.5-11



EZ Way is proud to be a US manufacturer of Safe Patient Handling products



PO Box 89, 807 E. Main • Clarinda, IA 51632 • 1-800-627-8940

Form #4-010
Rev. 08/18/2025